



POLYCARBONATE F SV Clear

Diam.	Sph.	Cyl.
65mm	+0.25 to +6.00	0.00 to -6.00
70mm	+0.25 to +6.00	0.00 to -2.00
70mm	0.00 to -6.00	0.00 to -6.00
70mm	-6.25 to -8.00	0.00 to -4.00
70mm	-8.25 to -10.00	0.00 to -2.00

POLYCARBONATE F SV PhotoPro Gray (Spin)

Diam.	Sph.	Cyl.
65mm	+0.25 to +4.00	0.00 to -4.00
70mm	+0.25 to +4.00	0.00 to -2.00
70mm	0.00 to -6.00	0.00 to -4.00

POLYCARBONATE F BIFOCAL (Clear & Photo Gray)

Diam.	Sph.	Add
Flat Top 70/28mm	+2.00 to -2.00	+1.00 to +3.50

POLYCARBONATE F PROGRESSIVE (Clear & Photo Gray)

Diam.	Sph.	Add
75/14mm	+2.00 to -2.00	+1.00 to +3.50
75/18mm	+2.00 to -2.00	+1.00 to +3.50



POLYCARBONATE SEMI-FINISHED SV

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Single vision	1.591	0.50, 1.50, 2.50, 3.50, 4.25, 5.25, 6.00, 7.50, 8.50	Clear	75	HC
Single vision	1.591	0.50, 1.50, 2.50, 3.50, 4.25, 5.25, 6.00, 7.50, 8.50	PhotoPro (Spin) Gray/Brown	75	HC
Single vision	1.591	0.50, 2.25, 4.25, 5.50, 7.50	Polarized Gray/Brown/G15	76	HC

POLYCARBONATE SEMI-FINISHED BIFOCAL

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Flat top	1.591	2.00, 4.00, 6.00	Clear	70/28	HC
Flat top	1.591	2.00, 4.00, 6.00	Gray	70/28	HC

POLYCARBONATE SEMI-FINISHED PROGRESSIVE

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Progressive	1.591	2.75, 4.75, 6.75	Clear	75	HC
Progressive	1.591	2.75, 4.75, 6.75	Gray	75	HC