

POLYCARBONATE UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|---------------|-------|--|------------|------------|---------|
| Single vision | 1.591 | 0.50, 1.50, 2.50, 3.50, 4.25, 5.25, 6.00, 7.50, 8.50 | Clear | 75 | HC |
| Single vision | 1.591 | 0.50, 1.50, 2.50, 3.50, 4.25, 5.25, 6.00, 7.50, 8.50 | Gray/Brown | 75 | HC |

POLYCARBONATE UV++ SEMI-FINISHED BIFOCAL

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|----------|-------|------------------|-------|------------|---------|
| Flat top | 1.591 | 2.00, 4.00, 6.00 | Clear | 70/28 | HC |
| Flat top | 1.591 | 2.00, 4.00, 6.00 | Gray | 70/28 | HC |
| Blended | 1.591 | 4.00, 6.00 | Clear | 70/28 | HC |

POLYCARBONATE UV++ SEMI-FINISHED PROGRESSIVE

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-------------|-------|------------------|-------|------------|---------|
| Progressive | 1.591 | 2.75, 4.75, 6.75 | Clear | 75 | HC |
| Progressive | 1.591 | 2.75, 4.75, 6.75 | Gray | 75 | HC |

CR-39 UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|--------------------------|-------|---|-----------|------------|---------|
| Single vision spherical | 1.499 | 0.00, 0.50, 1.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 70 | UC |
| Single vision spherical | 1.499 | 0.00, 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 75 | UC |
| Single vision aspherical | 1.499 | 0.00, 2.00, 3.00, 4.00, 6.00, 7.00, 8.00 | Clear | 75 | UC |
| Single vision spherical | 1.499 | 2.00, 4.025, 6.00, 7.00, 8.00 | Polarized | 77 | UC |
| Single vision spherical | 1.499 | 0.00, 1.00, 2.00, 3.00, 4.00, 5.25, 6.00, 7.25, 8.00, 9.00, 10.00 | Clear | 80 | UC |

CR-39 UV++ SEMI-FINISHED BIFOCAL

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-----------|-------|---|-------|------------|---------|
| Flat top | 1.499 | 0.50, 2.00, 4.00, 6.00, 8.00, 10.00 | Clear | 70/25 | UC |
| Flat top | 1.499 | 0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00, 10.00 | Clear | 70/28 | UC |
| Flat top | 1.499 | 0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00, 10.00 | Clear | 75/28 | UC |
| Flat top | 1.499 | 4.00, 6.00, 8.00 | Clear | 80/28 | UC |
| Round top | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | UC |
| Blended | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | UC |

CR-39 UV++ SEMI-FINISHED PROGRESSIVE

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-------------|-------|------------------------|-------|------------|---------|
| Progressive | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/12 | UC |
| Progressive | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/14 | UC |
| Progressive | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/17 | UC |
| Progressive | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/14 | UC |
| Progressive | 1.499 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/17 | UC |

1.56 UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|--------------------------|-------|---|------------------|------------|---------|
| Single vision spherical | 1.56 | 0.00, 0.50, 1.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 70 | HC |
| Single vision aspherical | 1.56 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 70 | HC |
| Single vision spherical | 1.56 | 0.00, 0.50, 1.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Photo gray/brown | 70 | HC |
| Single vision spherical | 1.56 | 0.00, 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 75 | HC |
| Single vision spherical | 1.56 | 0.00, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Photo gray/brown | 75 | HC |
| Single vision aspherical | 1.56 | 0.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00 | Clear | 75 | HC |
| Single vision spherical | 1.56 | 0.00, 1.00, 2.00, 3.00, 4.00, 5.25, 6.00, 7.25, 8.00, 9.00, 10.00 | Clear | 80 | HC |
| Single vision spherical | 1.56 | 0.00, 2.00, 4.00, 6.00, 8.00 | Photo gray/brown | 80 | HC |

1.56 UV++ SEMI-FINISHED BIFOCAL

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-----------|-------|-------------------------------------|------------------|------------|---------|
| Flat top | 1.56 | 0.50, 2.00, 4.00, 6.00, 8.00, 10.00 | Clear | 70/28 | HC |
| Flat top | 1.56 | 0.50, 2.00, 4.00, 6.00, 8.00, 10.00 | Photo gray/brown | 70/28 | HC |
| Flat top | 1.56 | 0.50, 2.00, 4.00, 6.00, 8.00, 10.00 | Clear | 75/28 | HC |
| Round top | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | HC |
| Round top | 1.56 | 2.00, 4.00, 5.00, 6.00, 8.00 | Photo gray/brown | 70/28 | HC |
| Blended | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | HC |
| Blended | 1.56 | 2.00, 4.00, 5.00, 6.00, 8.00 | Photo gray/brown | 70/28 | HC |

1.56 UV++ SEMI-FINISHED PROGRESSIVE

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-------------|-------|------------------------|------------------|------------|---------|
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/12 | HC |
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Photo gray/brown | 70/12 | HC |
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/14 | HC |
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Photo gray/brown | 70/14 | HC |
| Progressive | 1.56 | 2.00, 4.00, 5.00, 7.00 | Clear | 70/17 | HC |
| Progressive | 1.56 | 2.00, 4.00, 5.00, 7.00 | Photo gray/brown | 70/17 | HC |
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/14 | HC |
| Progressive | 1.56 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/17 | HC |

1.60 UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|--------------------------|-------|---|-----------|------------|---------|
| Single vision spherical | 1.60 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.25 | Clear | 70/75 | UC |
| Single vision spherical | 1.60 | 2.00, 4.00, 6.00, 8.00 | Polarized | 75 | UC/HC |
| Single vision aspherical | 1.60 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.25, 9.25, 10.25 | Clear | 70/75 | UC |
| Single vision spherical | 1.60 | 0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00 | Gray | 80 | UC/HC |

1.60 UV++ SEMI-FINISHED BIFOCAL

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-----------|-------|------------------------------|-------|------------|---------|
| Flat top | 1.60 | 0.50, 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | UC |
| Round top | 1.60 | 2.00, 4.00, 6.00, 8.00 | Clear | 70/28 | UC |

1.60 UV++ SEMI-FINISHED PROGRESSIVE

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-------------|-------|------------------------|-------|------------|---------|
| Progressive | 1.60 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/14 | UC |
| Progressive | 1.60 | 2.00, 4.00, 6.00, 8.00 | Clear | 75/17 | UC |

1.67 UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|--------------------------|-------|---|-----------|------------|---------|
| Single vision aspherical | 1.67 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.50 | Clear | 70 | UC |
| Single vision aspherical | 1.67 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.50 | Clear | 75 | UC |
| Single vision aspherical | 1.67 | 0.00, 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Grey | 75 | UC |
| Single vision aspherical | 1.67 | 2.00, 4.00, 6.00, 8.00 | Polarized | 75 | UC |
| Single vision aspherical | 1.67 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 80 | UC |

1.67 UV++ SEMI-FINISHED BIFOCAL

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|----------|-------|------------------------|-------|------------|---------|
| Flat top | 1.67 | 1.00, 3.00, 5.00, 7.00 | Clear | 75/28 | UC |

1.67 UV++ SEMI-FINISHED PROGRESSIVE

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|-------------|-------|------------------|-------|------------|---------|
| Progressive | 1.67 | 1.50, 3.50, 5.75 | Clear | 75/12 | UC/HC |

1.74 UV++ SEMI-FINISHED SV

| TYPE | INDEX | BASE CURVE | COLOR | BLANK SIZE | COATING |
|--------------------------|-------|---|-------|------------|---------|
| Single vision spherical | 1.74 | 0.50, 1.50, 2.25, 3.50, 5.00, 6.75, 8.50, 10.25 | Clear | 70/75 | UC |
| Single vision aspherical | 1.74 | 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00 | Clear | 70/75 | UC |