

POLYCARBONATE UV++ F SV (Clear & Photo Gray)

Diam.	Sph.	Cyl.
65mm	+0.25 to +6.00	0.00 to -6.00
70mm	+0.25 to +6.00	0.00 to -2.00
70mm	+0.00 to -6.00	0.00 to -6.00
70mm	-6.25 to -8.00	0.00 to -4.00
70mm	-8.25 to -10.00	0.00 to -2.00

POLYCARBONATE UV++ F BIFOCAL (Clear & Photo Gray)

Bifocal	Diam.	Sph.	Add
Flat Top	70/28mm	+2.00 to -2.00	+1.00 to +3.50
Blended	70/28mm	+2.00 to -2.00	+1.00 to +3.50

POLYCARBONATE UV++ F PROGRESSIVE (Clear & Photo Gray)

Progressive	Diam.	Sph.	Add
Short corridor	75/14mm	+2.00 to -2.00	+1.00 to +3.50
Long corridor	75/18mm	+2.00 to -2.00	+1.00 to +3.50

CR-39 UV++ F SV

Diam.	Sph.	Cyl.
55mm	+0.25 to +8.00	0.00 to -2.00
60mm	+0.25 to +8.00	0.00 to -4.00
60mm	+8.25 to +12.00	0.00 to -2.00
65mm	+6.00 to -6.00	0.00 to -2.00
65mm	+6.00 to -6.00	0.00 to -4.00
70mm	+6.00 to -20.00	0.00 to -2.00
70mm	-0.00 to -8.00	-2.25 to -4.00
70mm	-0.00 to -6.00	-4.25 to -6.00

CR-39 UV++ F BIFOCAL

Bifocal	Diam.	Sph.	Add
Flat Top	70/28mm	+2.00 to -2.00	+1.00 to +3.50
Round Top	70/24mm	+2.00 to -2.00	+1.00 to +3.50
Round Top	70/28mm	+2.00 to -2.00	+1.00 to +3.50
Blended	70/24mm	+2.00 to -2.00	+1.00 to +3.50
Blended	70/28mm	+2.00 to -2.00	+1.00 to +3.50

CR-39 UV++ F PROGRESSIVE

Progressive	Diam.	Sph.	Add
Short corridor	70/12mm	+2.00 to -2.00	+1.00 to +3.50
Short corridor	70/14mm	+2.00 to -2.00	+1.00 to +3.50
Short corridor	75/14mm	+2.00 to -2.00	+1.00 to +3.50
Long corridor	75/17mm	+2.00 to -2.00	+1.00 to +3.50

CR-39 UV++ PLANO LENS

Diam.	Thickness	Base Curve
71mm	1.8mm	4B/5B/6B
71mm	2.0mm	4B/5B/6B
73mm	1.8mm	4.5/6B
73mm	2.0mm	4.5/6B
75mm	1.8mm	0/2/4/5/6B
75mm	2.0mm	0/2/4/5/6B
80mm	1.8mm	6B
80mm	2.0mm	6B

1.56 UV++ F SV

Diam.	Sph.	Cyl.
65mm	+0.25 to +8.00	0.00 to -2.00
65mm	+0.25 to +8.00	0.00 to -4.00
72mm	-6.00 to -20.00	0.00 to -2.00
72mm	-0.00 to -8.00	-2.25 to -4.00
72mm	-0.00 to -6.00	-4.25 to -6.00

1.56 UV++ F BIFOCAL

Bifocal	Diam.	Sph.	Add
Flat Top	70/28mm	+2.00 to -2.00	+1.00 to +3.50
Round Top	70/24mm	+2.00 to -2.00	+1.00 to +3.50
Round Top	70/28mm	+2.00 to -2.00	+1.00 to +3.50
Blended	70/24mm	+2.00 to -2.00	+1.00 to +3.50
Blended	70/28mm	+2.00 to -2.00	+1.00 to +3.50

1.56 UV++ F PROGRESSIVE

Progressive	Diam.	Sph.	Add
Short corridor	70/12mm	+2.00 to -2.00	+1.00 to +3.50
Short corridor	70/14mm	+2.00 to -2.00	+1.00 to +3.50
Long corridor	70/17mm	+2.00 to -2.00	+1.00 to +3.50

1.60 UV++ F SV

Diam.	Sph.	Cyl.
65mm	+0.25 to +6.00	0.00 to -2.00
70mm	-0.00 to -8.00	0.00 to -4.00
70mm	-8.25 to -18.00	0.00 to -2.00
75mm	-0.00 to -6.00	0.00 to -2.00

1.60 UV++ PHOTOCROMIC F SV

Diam.	Sph.	Cyl.
65mm	+0.25 to +6.00	0.00 to -2.00
70mm	-0.00 to -8.00	0.00 to -4.00
70mm	-8.25 to -18.00	0.00 to -2.00
75mm	-0.00 to -6.00	0.00 to -2.00

1.67 UV++ F SV

Diam.	Sph.	Cyl.
65mm	+0.25 to +6.00	0.00 to -2.00
70mm	-0.00 to -12.00	0.00 to -2.00
70mm	-12.50 to -15.00	0
70mm	-0.00 to -12.00	0.00 to -4.00
75mm	-0.00 to -6.00	0.00 to -2.00

Diam.	Sph.	Cyl.
70mm	-6.25 to -10.00	0.00 to -4.00
70mm	-10.25 to -12.00	0.00 to -2.00
70mm	-12.50 to -15.00	0.00
75mm	-0.00 to -6.00	0.00 to -4.00